

Silhouette™ Arm Undersleeve Sizing Chart - Standard Size

Facility Name / Phone Number

Contact / Fax Number

Billing and/or Shipping Address Line 1

Patient's Name

Billing and/or Shipping Address Line 2

Purchase Order No. or Reference / Account No.

City, State, Zip

Shipping Method					
<input type="checkbox"/> Ground	<input type="checkbox"/> 3 Day	<input type="checkbox"/> 2 Day	<input type="checkbox"/> Next Day	<input type="checkbox"/> () Early am	<input type="checkbox"/> () 10:30

Please take measurements in centimeters

1. Take circumference measurements when limb is at its smallest with patients arm slightly bent at 35 degrees.
2. Measure arm circumference according to chart below, starting at palm of hand up to axilla.
3. Medial Length from wrist to axilla is: _____ (Standard sizes fit 39cm - 46cm)
4. Medial length from wrist to elbow crease is: _____ (Standard sizes fit 19cm - 24cm)
5. Be sure measurements are accurate. When on borderline between sizes, consider amount of edema present and anticipated reduction. If significant edema is present and reduction is anticipated, select smaller size. If edema is minimal, select larger size. Also consider desired fit. Measurements at maximum of size range will return a sleeve with a snug fit, while measurements at minimum of size range will yield a roomier fitting sleeve.
6. If measurements fall outside these parameters, measure for a custom Silhouette Undersleeve (separate form).

*A reinforced seam can be cut open to effectively shorten sleeve's fit by 2.5cm. See instructions for use.

Circumference Measurements			Small	Medium	Large
Position	Left	Right	centimeters	centimeters	centimeters
Palm of Hand			16 - 21	17 - 23	19 - 24
Wrist			12 - 23	15 - 26	21 - 32
5cm from wrist			13 - 24	17 - 28	23 - 34
15cm from wrist			18 - 26	23 - 32	29 - 38
25cm from wrist			22 - 28	28 - 35	34 - 41
35cm from wrist			23 - 31	30 - 38	36 - 44

